


# Tap Density Tester



## USP 1 and 2 (EDT-1020)



2  
Methods

USP Method 1  
and USP Method 2



Simultaneous Rotating  
and Tapping Motion



User Friendly  
Operation



LAN (optional)



Snaplock™  
Mechanism



Free Drop of  
the Cylinder



Print




Report

### Key Features

- Supports USP 1 (300 taps per min) and USP 2 (250 taps per min) and ASTM test methods
- Simultaneous rotating and tapping motion assures an evenly packed surface
- Calculation of test results like tapped density, compressibility index and Hausner ratio
- Printer port for documenting test results as per GMP/GLP standards
- LAN connectivity for data transfer to a computer

USP 1 (ETD-10) & USP 2 (ETD-20)  
is also available




Tap Density Test Report

# Tap Density Tester



## Optional Accessories



Acoustic cabinet

This reduces the sound level to 71 dB to meet laboratory standards



Adapter for 50 ml cylinder



Adapter for 10 ml cylinder