

Automatic Manipulator

for Cell Factory, CellSTACK, HYPERStack (3 mounting type)

WAS-111

Appearance




Features

- With this equipment, **stirring and liquid handling** can be performed simultaneously in multiple stackable cell culture vessels.*
- Up to this equipment can mount **3 pieces of 40(or 52) layers stackable cell culture vessels**.
 [Options can be added :
 Mounting **6** stackable cell culture vessels (10 layer Cell Factory , 10 layer CellSTACK ,
 or 36 layer HYPERStack) at maximum.]
- **Vertical, horizontal, diagonal action** in addition to complicated 3D action are possible.

*Cell Factory, CellStack , HYPERStack

Multi-Language HMI



•Language selectable.

